



Sulfonamides Residue Rapid Test Strip (Milk) (DTS008)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Size	40T
Intended Use	Sulfonamides device is for rapid test to qualitatively detect the Sulfonamides in milk sample at the sensitivity of 50-200 µg/kg. It takes approx. 5~8 min.
Storage	Store at 15-25°C, DO NOT FREEZE or use beyond the expiration date. The shelf life is 12 months.
Sensitivity	<p>SMZ.....50 ppb SD.....100 ppb SMM.....100 ppb ST..... 100 ppb SM1..... 200 ppb</p> <p>SMZ: sulfamethizole SD: sulfadiazine SMM: Sulfamonomethoxine ST: Sulfathiazole SM1: sulfamerazine</p>