



Salbutamol Residue Rapid Test Strip (Meat) (DTS010)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Size 40T

Intended Use

The Salbutamol Residue Rapid Inspection Device is a rapid, one step test for the qualitative detection of Salbutamol in meat samples. The whole time of assay is only approx. 5~8 min.

Storage

Store at 15-25°C, DO NOT FREEZE or use beyond the expiration date. The shelf life is 12 months.

Sensitivity

To acquire the exact sensitivity, reduplicative experiment has been done on the sample containing 5 µg/kg Salbutamol.

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 Fax: 1-631-938-8221

Email:info@creative-diagnostics.com