



Tetracycline Residue Rapid Test Strip (Milk) (DTS019)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Size	40T
Intended Use	Tetracycline device is for rapid test to qualitatively detect the Tetracycline in milk sample at the sensitivity of 10 µg/kg. It only takes approx. 5~10 min.
Storage	Store at 4-30°C, DO NOT FREEZE or use beyond the expiration date. The shelf life is 12 months.
Sensitivity	To acquire the exact sensitivity, reduplicative experiment has been done. Tetracycline.....10ppb Oxybicycline.....10ppb Deoxycycline.....15ppb Chlorotetracycline.....30 ppb