



Tetracycline Residue Rapid Test Strip (Tissue) (DTS020)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Size	40T
Intended Use	The Tetracycline Residue Rapid Test Device is a rapid, one step test for the qualitative (or quantitative) detection of most common Tetracyclines residues in tissues at the sensitivity of 10 µg/kg (10 ppb). The total assay time takes approx. 10min.
Storage	Store at 4-30°C, DO NOT FREEZE or use beyond the expiration date. The shelf life is 12 months.