



Penicilin Residue Rapid Test Strip (Honey) (DTS033)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Size	40T
Intended Use	Penicilin Residue Rapid Test Device is a rapid, one step test for the qualitative detection of Penicilin in honey samples. It takes only approx. 5~10 min.
Storage	Store at 4-30°C, DO NOT FREEZE or use beyond the expiration date. The shelf life is 12 months.
Sensitivity	To acquire the exact sensitivity, reduplicative experiment has been done on the sample containing 4 µg/kg Penicilin.
Specificity	No cross-reaction with Eepharmycin, Chloramphicol, Streptomycin, or sulfametazine.