



## Pseudotyped GFP HPV35 (PSVG-HPV35)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	HPV35 Pseudovirus are used to test the ability of serum and antibodies to neutralize the infectivity of HPV35. HPV35 Pseudovirus are produced in 293FT cells using three separate plasmids, encoding the HPV L1 protein, the HPV L2 protein, and a reporter gene. HPV35 Pseudovirus display antigenically correct L1 and L2 protein pseudotyped on replication-incompetent virus particles. Pseudovirus are capable of a single round of infection and carry a plasmid that expresses GFP optical reporter gene upon infection.
<b>Applications</b>	Useful for studying papillomavirus assembly, entry, and neutralization. The recommended amount of pseudovirus per well can be found in the Certificate of Analysis for each lot of pseudovirus.
<b>Size</b>	100 µl
<b>Storage</b>	Store at -80°C. Multiple freeze/thaw cycles not recommended. When using the virus, transfer the virus from the -80 ° C refrigerator and melt it in an ice bath.
<b>Ship</b>	Frozen on dry ice