



Anti-CD4 monoclonal antibody, clone YTS177.9 [R-PE] (CABT-45116RM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview Rat anti Mouse CD4 antibody, clone YTS177.9 reacts with the Mouse CD4 antigen, non polymorphic epitope. Rat anti Mouse CD4 antibody, Rat anti Mouse CD4 antibody, clone YTS177.9 reacts with CD4 transfectants. Rat anti Mouse CD4 antibody, clone YTS177.9 also blocks MHC-II dependant T-cell responses in vitro and in vivo and induces tolerance. Flow Cytometry Use 10ul of the suggested working dilution to label 10⁶ cells or 100ul of whole blood. The Fc region of monoclonal antibodies may bind non-specifically to cells expressing low affinity Fc receptors.

Specificity	CD4
Immunogen	Mouse spleen cells
Isotype	IgG2a
Source/Host	Rat
Species Reactivity	Mouse
Clone	YTS177.9
Conjugate	PE
Applications	FC
Format	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilised
Size	100 tests
Preservative	0.09% Sodium Azide
Storage	Prior to reconstitution store at +4°C. Following reconstitution store at +4°C. DO NOT FREEZE.

This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Cd4 CD4 antigen [Mus musculus (house mouse)]
Official Symbol	CD4
Synonyms	CD4; CD4 antigen; L3T4; Ly-4; T-cell surface glycoprotein CD4; T-cell surface antigen T4/Leu-3; T-cell differentiation antigen L3T4; T-cell surface glycoprotein CD4 precursor (T-cell surface antigen T4/Leu-3) (T-cell differentiation antigen L3T4);
Entrez Gene ID	12504
Protein Refseq	NP_038516
UniProt ID	P06332
Chromosome Location	6 F2; 6 59.17 cM
Pathway	Adaptive Immune System; Alpha-defensins; Antigen processing and presentation; Cell adhesion molecules (CAMs); Costimulation by the CD28 family; Cytokines and Inflammatory Response (BioCarta); Defensins; Downstream TCR signaling;
Function	coreceptor activity; enzyme binding; glycoprotein binding; protein binding; protein homodimerization activity; protein kinase binding; zinc ion binding;