



Anti-CD8A monoclonal antibody, clone 53-6.7 [R-PE/Cy5®] (CABT-45306RM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

Rat anti Mouse CD8 alpha antibody, clone 53-6.7 detects mouse CD8a, the alpha chain of the CD8 antigen. CD8 is expressed on T-cells, and exists either a homodimer composed of two alpha chains, or a heterodimer composed of one alpha and one beta chain. CD8 functions as a co-receptor with the T-cell receptor and mediates efficient cell to cell interactions within the immune system. The CD8 alpha chain specifically binds to class I MHC molecules. Rat anti Mouse CD8 alpha antibody, clone 53-6.7 has been reported to block antigen presentation via MHC class I and inhibit T cell responses to IL-2. This antibody has also been used for depletion of CD8a+ cells. Flow Cytometry Use 5ul of the suggested working dilution to label 1x10⁶ cells in 100ul. The Fc region of monoclonal antibodies may bind non-specifically to cells expressing low affinity Fc receptors.

Specificity	CD8A
Immunogen	Mouse thymus or spleen.
Isotype	IgG2a
Source/Host	Rat
Species Reactivity	Mouse
Clone	53-6.7
Conjugate	PE/Cy5
Applications	FC
Format	Purified IgG conjugated to R. Phycoerythrin (RPE) Cy5 - liquid
Size	500 µl

Preservative	0.09% Sodium Azide
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Cd8a CD8 antigen, alpha chain [Mus musculus (house mouse)]
Official Symbol	CD8A
Synonyms	CD8A; CD8 antigen, alpha chain; Ly-2; Ly-B; Ly-35; Lyt-2; BB154331; T-cell surface glycoprotein CD8 alpha chain; T-cell surface glycoprotein Lyt-2; Lyt-2.1 lymphocyte differentiation antigen (AA at 100);
Entrez Gene ID	12525
Protein Refseq	NP_001074579
UniProt ID	P01731
Chromosome Location	6 C; 6 32.14 cM
Pathway	Adaptive Immune System; Antigen processing and presentation; Cell adhesion molecules (CAMs); Hematopoietic cell lineage; Immune System; Immunoregulatory interactions between a Lymphoid and a non-Lymphoid cell; Primary immunodeficiency; T Cell Receptor Signaling Pathway;
Function	protein binding; protein homodimerization activity; protein kinase binding;