



Anti-Human IgG monoclonal antibody, clone A57H (CABT-48847MH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview Mouse anti Human IgG (CH3 Domain) antibody, clone A57H reacts with the CH3 domain of human IgG Immunohistology This product does not require protein digestion pre-treatment of paraffin sections. This product does not require antigen retrieval using heat treatment prior to staining of paraffin sections.

Specificity	IgG
Immunogen	pFc' from purified protein of Heavy Chain Disease.
Isotype	IgM
Source/Host	Mouse
Species Reactivity	Human
Clone	A57H
Conjugate	Unconjugated
Applications	ELISA, HA, WB
Format	Purified IgM - liquid
Size	200 µg
Preservative	See individual product datasheet
Storage	Store undiluted at -20°C

GENE INFORMATION

Gene Name	IGHG1 immunoglobulin heavy constant gamma 1 (G1m marker) [Homo sapiens (human)]
Official Symbol	IGHG1
Synonyms	IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IGG;
Entrez Gene ID	3500
UniProt ID	P01860
Chromosome Location	14q32.33
Pathway	IL4-mediated signaling events;
Function	antigen binding; protein binding;