



Anti-SMAD6 polyclonal antibody (CABT-BL001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein of human SMAD6
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB
Molecular Weight	53kDa
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	The protein encoded by this gene belongs to the SMAD family of proteins, which are related to Drosophila 'mothers against decapentaplegic' (Mad) and C. elegans Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways. This protein functions in the negative regulation of BMP and TGF-beta/activin-signalling.
--------------	---

Multiple transcript variants encoding different isoforms have been found for this gene.

GENE INFORMATION

Entrez Gene ID [4091](#)

UniProt ID [O43541](#)
