



Anti-STX8 polyclonal antibody (CABT-BL003)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Immunogen | Recombinant protein of human STX8 |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human, Mouse |
| Purification | Affinity purification |
| Conjugate | Unconjugated |
| Applications | WB |
| Molecular Weight | 27kDa |
| Format | Liquid |
| Size | 100 µl |
| Buffer | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3. |
| Preservative | 0.02% Sodium Azide |
| Storage | Store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|---------------------|--|
| Introduction | The gene is a member of the syntaxin family. The encoded protein is involved in protein trafficking from early to late endosomes via vesicle fusion and exocytosis. A related pseudogene has been identified on chromosome 12. Alternative splicing results in multiple transcript variants. |
|---------------------|--|

GENE INFORMATION

Entrez Gene ID [9482](#)

UniProt ID [Q9UNK0](#)
