



Anti-DLGAP5 polyclonal antibody (CABT-BL012)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein of human DLGAP5
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF
Molecular Weight	95kDa
Format	Liquid
Size	100 µl
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	DLG7 is a potential cell cycle regulator that may play a role in carcinogenesis of cancer cells and mitotic phosphoprotein regulated by the ubiquitin-proteasome pathway. The protein is also
--------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

a key regulator of adherens junction integrity and differentiation that may be involved in CDH1-mediated adhesion and signaling in epithelial cells.

GENE INFORMATION

Entrez Gene ID	9787
----------------	----------------------

UniProt ID	Q15398
------------	------------------------
