



Rabbit Anti-Histone H2B monoclonal antibody, clone TE31-74 (CABT-BL8535)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Histone H2B
Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	TE31-74
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	WB, ICC/IF, IHC, IP, FC
Molecular Weight	14 kDa
Cellular Localization	Nucleus, Chromosome.
Positive Control	NIH/3T3, PC12, Hela, A431, human tonsil tissue, human liver tissue, human breast carcinoma tissue, mouse liver tissue, mouse testis tissue, mouse colon tissue.
Format	Liquid
Size	100 µl
Buffer	1×TBS (pH7.4), 1% BSA, 40% Glycerol.

