



Anti-beta Galactosidase (LacZ) polyclonal antibody (CABT-BL8691)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with beta-galactosidase as determined by IEP and ELISA techniques.
Target	beta Galactosidase
Isotype	IgY
Source/Host	Chicken
Species Reactivity	N/A
Purification	Immunoaffinity purified
Conjugate	Unconjugated
Applications	ELISA, IEP, IHC, WB This antibody reacts with beta-galactosidase as determined by IEP and ELISA techniques. It is suitable for blotting, ELISA and IHC applications. Optimal working dilutions should be determined experimentally by the investigator.
Format	Liquid
Size	100 µg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 2 years