



# Anti-Vibrio Cholera Toxin B monoclonal antibody (CABT-BL8728)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Cholera Toxin
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Cholera
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Lyophilized
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	None
Storage	4°C for 1 year