



Anti-CSFV Structural Glycoprotein monoclonal antibody, clone 1D2 (CABT-BL8758)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	CSFV-E2
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	CSFV
Clone	1D2
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Lyophilized
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	None
Storage	4°C for 1 year