



Anti-Dengue Envelope Protein monoclonal antibody (CABT-BL8778)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Dengue Envelope
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Dengue
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IF, WB
Format	Liquid
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 1 year