



Anti-Giardia lamblia monoclonal antibody, clone 2G8 (CABT-BL8822)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Giardia-CWP
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Giardia
Clone	2G8
Purification	Purified Protein
Conjugate	Unconjugated
Applications	ELISA, IF, WB
Format	Lyophilized
Size	1 mg
Buffer	PBS
Preservative	0.1% Sodium Azide
Storage	-20°C for 1 year