



Anti-Haptoglobin monoclonal antibody (CABT-BL8836)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Haptoglobin
Source/Host	Mouse
Species Reactivity	Horse
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IP, WB
Format	Lyophilized
Size	1 mg
Buffer	Phospate buffered saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	4°C for 1 year