



Rabbit Anti-PAK1, PAK2, PAK3 monoclonal antibody, clone TV14-81 (CABT-L654)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Phospho-PAK1(S144)+PAK2(S141)+PAK3(S139)
Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	TV14-81
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	WB, ICC/IF, IHC, IP, FC
Molecular Weight	65 kDa
Cellular Localization	Cytoplasm, Cell membrane.
Positive Control	NIH/3T3, SH-SY-5Y, HeLa, human liver cancer tissue, mouse brain tissue, human colon cancer tissue.
Format	Liquid
Size	100 µl
Buffer	1×TBS (pH7.4), 1% BSA, 40% Glycerol.

Preservative 0.05% Sodium Azide

Storage Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Cytoplasm, Cell membrane.
