



Anti-C1QC polyclonal antibody (CPBT-65027SH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep anti Human C1q antibody recognizes human C1q, a sub-component of the C1 complement component.
Specificity	C1QC
Immunogen	Human C1q, purified from plasma
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Conjugate	FITC
Applications	IHC-Fr
Format	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Size	1 ml
Preservative	0.09% Sodium Azide
Storage	in frost free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name [C1QC complement component 1, q subcomponent, C chain \[Homo sapiens \(human\) \]](#)

Official Symbol	C1QC
Synonyms	C1QC; complement component 1, q subcomponent, C chain; C1QG; C1Q-C; complement C1q subcomponent subunit C; complement component 1, q subcomponent, gamma polypeptide;
Entrez Gene ID	714
Protein Refseq	NP_001107573
UniProt ID	P02746
Chromosome Location	1p36.11
Pathway	Chagas disease (American trypanosomiasis); Classical antibody-mediated complement activation; Complement Activation, Classical Pathway; Complement and Coagulation Cascades; Complement and coagulation cascades; Complement cascade; Creation of C4 and C2 activators; Immune System;