



Anti-C1QC polyclonal antibody (CPBT-65028SH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep anti Human C1q antibody recognizes human C1q, a sub-component of the C1 complement component.
Specificity	C1QC
Immunogen	Human C1q, purified from plasma
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA; FC; ID
Format	Purified IgG - liquid
Size	1 ml
Preservative	0.09% Sodium Azide
Storage	in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name [C1QC complement component 1, q subcomponent, C chain \[Homo sapiens \(human\) \]](#)

Official Symbol	C1QC
Synonyms	C1QC; complement component 1, q subcomponent, C chain; C1QG; C1Q-C; complement C1q subcomponent subunit C; complement component 1, q subcomponent, gamma polypeptide;
Entrez Gene ID	714
Protein Refseq	NP_001107573
UniProt ID	P02746
Chromosome Location	1p36.11
Pathway	Chagas disease (American trypanosomiasis); Classical antibody-mediated complement activation; Complement Activation, Classical Pathway; Complement and Coagulation Cascades; Complement and coagulation cascades; Complement cascade; Creation of C4 and C2 activators; Immune System;