



Anti-Human IgG Fc polyclonal antibody (CPBT-65491GH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This product is specific for human IgG (Fc). Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-goat serum, human IgG, human IgG F(c)and human serum.
Specificity	IgG
Immunogen	Human IgG (Fc) fragment
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA; IP; WB
Format	Purified IgG - liquid
Size	2 mg
Preservative	0.09% Sodium Azide
Storage	in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name [IGHG1 immunoglobulin heavy constant gamma 1 \(G1m marker\) \[Homo sapiens \(human\) \]](#)

Official Symbol	IGHG1
Synonyms	IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IGG;
Entrez Gene ID	3500
UniProt ID	P01860
Chromosome Location	14q32.33
Pathway	IL4-mediated signaling events;
Function	antigen binding; protein binding;