



Anti-Human IgG Fc polyclonal antibody (CPBT-65512GH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This product recognizes human IgG (Fc). This reagent is available in bulk quantities and may be made available in other formats.
Specificity	IgG
Immunogen	Purified human IgG (Fc)
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Conjugate	Unconjugated
Applications	Turbidimetry
Format	Serum - liquid
Size	1 ml
Preservative	0.09% Sodium Azide
Storage	in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name [IGHG1 immunoglobulin heavy constant gamma 1 \(G1m marker\) \[Homo sapiens \(human\) \]](#)

Official Symbol	IGHG1
Synonyms	IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IGG;
Entrez Gene ID	3500
UniProt ID	P01857
Chromosome Location	14q32.33
Pathway	IL4-mediated signaling events;
Function	antigen binding; protein binding;