



Anti-Human IgG polyclonal antibody [HRP] (CPBT-65531GH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	IgG
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Conjugate	HRP
Applications	IHC-Fr; ELISA
Preservative	None

GENE INFORMATION

Gene Name	IGKC immunoglobulin kappa constant [Homo sapiens (human)]
Official Symbol	IGKC
Synonyms	IGKC; immunoglobulin kappa constant; Km; HCAK1; IGKCD; IGG;
Entrez Gene ID	3514
UniProt ID	P01834
Chromosome Location	2p12
Pathway	Adaptive Immune System; Antigen activates B Cell Receptor (BCR) leading to generation of second messengers; Binding and Uptake of Ligands by Scavenger Receptors; Classical

antibody-mediated complement activation; Complement cascade; Creation of C4 and C2 activators; FCERI mediated Ca+2 mobilization; FCERI mediated MAPK activation;

Function

antigen binding;
