



Anti-C3 polyclonal antibody (CPBT-65830GH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This product recognises human C3, a secreted protein that plays a central role in the activation of the complement system. This reagent is available in bulk quantities and may be made available in other formats. Please enquire for further details.
Specificity	C3
Immunogen	Purified human C3
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Conjugate	Unconjugated
Applications	Turbidimetry
Format	Serum - liquid
Size	1 ml
Preservative	0.09% Sodium Azide
Storage	in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name [C3 complement component 3 \[Homo sapiens \(human\)\]](#)

Official Symbol	C3
Synonyms	C3; complement component 3; ASP; C3a; C3b; AHUS5; ARMD9; CPAMD1; HEL-S-62p; complement C3; prepro-C3; C3a anaphylatoxin; complement component C3; complement component C3a; complement component C3b; acylation-stimulating protein cleavage product; epididymi
Entrez Gene ID	718
Protein Refseq	NP_000055
UniProt ID	P01024
Chromosome Location	19p13.3-p13.2
Pathway	Activation of C3 and C5; Adaptive Immune System; Alternative complement activation; Chagas disease (American trypanosomiasis); Class A/1 (Rhodopsin-like receptors); Complement Activation, Classical Pathway; Complement and Coagulation Cascades; Complement and coagulation cascades;
Function	C5L2 anaphylatoxin chemotactic receptor binding; endopeptidase inhibitor activity; protein binding; receptor binding;