



Anti-Rat IgE polyclonal antibody [HRP] (CPBT-65987SR)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	IgE
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Rat
Conjugate	HRP
Applications	IHC-Fr; ELISA; WB
Size	500 µg
Preservative	None

GENE INFORMATION

Gene Name	Ighe Immunoglobulin heavy constant epsilon [Mus musculus (house mouse)]
Official Symbol	Ighe
Synonyms	IGHE; Immunoglobulin heavy constant epsilon; Gm900; IGE;
Entrez Gene ID	380792
UniProt ID	P06336
Chromosome Location	12; 12 F1

Functionmolecular_function
