



Anti-Albumin polyclonal antibody [HRP] (CPBT-66071SR)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep anti Rat albumin antibody recognizes rat albumin, a soluble protein which constitutes approximately 50% of blood serum protein, acting as a regulator of blood colloidal osmotic pressure and hence blood volume. albumin is also an important transporter of substances within the bloodstream, such as fatty acids, thyroid hormones, metal ions and steroids.
Specificity	ALB
Target	Albumin
Immunogen	Native albumin from rat serum
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Rat, Bovine
Conjugate	HRP
Applications	ELISA
Format	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Size	1 ml
Preservative	See individual product datasheet
Storage	in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Alb albumin [Rattus norvegicus (Norway rat)]
Official Symbol	ALB
Synonyms	ALB; albumin; Alb1; Albza; serum albumin;
Entrez Gene ID	24186
Protein Refseq	NP_599153
UniProt ID	P02770
Chromosome Location	14p22
Pathway	Bile acid and bile salt metabolism; Binding and Uptake of Ligands by Scavenger Receptors; HDL-mediated lipid transport; Hemostasis; Lipid digestion, mobilization, and transport; Lipoprotein metabolism; Metabolism; Metabolism of lipids and lipoproteins;
Function	DNA binding; chaperone binding; drug binding; enzyme binding; fatty acid binding; contributes_to oxygen binding; pyridoxal phosphate binding; toxic substance binding; zinc ion binding;