



Anti-ALPI polyclonal antibody (CPBT-67302GC)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This product recognises calf intestinal alkaline phosphatase.
Specificity	ALPI
Immunogen	Highly pure calf intestinal alkaline phosphatase
Isotype	IgG
Source/Host	Goat
Species Reactivity	Bovine
Conjugate	Unconjugated
Applications	IHC-Fr; ELISA; IHC-P
Format	Serum - liquid
Size	1 ml
Preservative	0.09% Sodium Azide
Storage	in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	ALPI alkaline phosphatase, intestinal [Bos taurus (cattle)]
Official Symbol	ALPI

Synonyms	ALPI; IAP; intestinal-type alkaline phosphatase; intestinal alkaline phosphatase V;
Entrez Gene ID	280993
Protein Refseq	NP_776412
UniProt ID	P19111
Chromosome Location	2q45
Pathway	Folate biosynthesis; Metabolic pathways.
Function	alkaline phosphatase activity; magnesium ion binding; zinc ion binding;