



# Anti-Albumin polyclonal antibody (CPBT-67345DS)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Donkey anti Sheep albumin antibody recognizes sheep albumin, a major plasma protein that functions in the regulation of colloidal osmotic pressure.
<b>Specificity</b>	ALB
<b>Target</b>	Albumin
<b>Immunogen</b>	Sheep albumin, essentially globulin free (approximately 99%) prepared from fraction V powder.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Donkey
<b>Species Reactivity</b>	Sheep
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Serum - liquid
<b>Size</b>	5 ml
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">ALB albumin [ Ovis aries (sheep) ]</a>
<b>Official Symbol</b>	ALB
<b>Synonyms</b>	ALB; serum albumin; pre-pro serum albumin;
<b>Entrez Gene ID</b>	<a href="#">443393</a>
<b>Protein Refseq</b>	<a href="#">NP_001009376</a>
<b>UniProt ID</b>	P14639
<b>Chromosome Location</b>	chromosome: 6