



## Anti-HEY1 polyclonal antibody (CPBT-67678RH)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Specificity	HEY1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-P; WB
Size	25 µg
Preservative	None

### GENE INFORMATION

Gene Name	<a href="#">HEY1 hes-related family bHLH transcription factor with YRPW motif 1 [ Homo sapiens (human) ]</a>
Official Symbol	HEY1
Synonyms	HEY1; hes-related family bHLH transcription factor with YRPW motif 1; CHF2; OAF1; HERP2; HESR1; HRT-1; hHRT1; BHLHb31; hairy/enhancer-of-split related with YRPW motif protein 1; HES-related repressor protein 1; HES-related repressor protein 2; basic helix
Entrez Gene ID	<a href="#">23462</a>
Protein Refseq	<a href="#">NP_001035798.1</a>

<b>UniProt ID</b>	Q9Y5J3
<b>Chromosome Location</b>	8q21
<b>Pathway</b>	Constitutive Signaling by NOTCH1 HD+PEST Domain Mutants; Constitutive Signaling by NOTCH1 PEST Domain Mutants; Delta-Notch Signaling Pathway; Disease; FBXW7 Mutants and NOTCH1 in Cancer;
<b>Function</b>	RNA polymerase II core promoter sequence-specific DNA binding transcription factor activity; microsatellite binding; protein binding; protein binding transcription factor activity