



Rabbit Anti-Mouse IgG1 polyclonal antibody [HRP] (CPBT-68016RM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Mouse IgG1 heavy chain
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Affinity purification
Conjugate	HRP
Applications	WB, ELISA, IHC Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Reconstitution	Restore with deionized water (or equivalent)
Format	Lyophilized
Size	1 mg
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 0.01% Gentamicin Sulfate.
Storage	Store at 4°C prior to restoration For extended storage aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing.

GENE INFORMATION

Official Symbol	IGHG1
Synonyms	IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IgG1; Igh-4; VH7183;
Entrez Gene ID	16017
UniProt ID	P01868