



Anti-Mouse IgG polyclonal antibody [R-PE] (CPBT-68109GA)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	OBT031R is specific for mouse IgG heavy and light chains. Minimal cross reactivity has been observed with bovine, horse, human, rabbit and pig serum proteins.
Specificity	IgG
Immunogen	Mouse IgG whole molecule
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	PE
Applications	FC; IF
Size	1 ml
Preservative	0.01% Sodium Azide
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Ighg2b immunoglobulin heavy constant gamma 2B [Mus musculus (house mouse)]
Official Symbol	IGHG2B

Synonyms	IGHG2B; immunoglobulin heavy constant gamma 2B; IgG2b; Igh-3; gamma2b; IGG;
Entrez Gene ID	16016
UniProt ID	P01845
Chromosome Location	12 F2; 12 62.1 cM
Function	antigen binding;