



Anti-Mouse IgG polyclonal antibody [AP] (CPBT-68110GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview This product recognises Mouse IgG. Cross-reactivity with Bovine, Chicken, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat or Sheep immunoglobulin has been minimised by solid phase absorption.

Specificity IgG

Immunogen Mouse IgG whole molecule

Isotype IgG

Source/Host Goat

Species Reactivity Mouse

Conjugate AP

Applications IHC-Fr; ELISA; WB

Format Purified IgG conjugated to Alkaline Phosphatase - liquid

Size 1 mg

Preservative 0.1% Sodium Azide

Storage Store at +4°C. DO NOT FREEZE. This product should be stored undiluted.

GENE INFORMATION

Gene Name [Ighg1 immunoglobulin heavy constant gamma 1 \(G1m marker\) \[Mus musculus \(house mouse\) \]](#)

Official Symbol	IGHG1
Synonyms	IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IgG1; Igh-4; VH7183; IGG;
Entrez Gene ID	16017
UniProt ID	P01837
Chromosome Location	12 F1; 12 62.1 cM
Function	antigen binding;