



S. hemolyticus Streptolysin O (DAG-T2696)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Purified native Streptolysin O protein for use in IVD diagnostics
Species	Streptococcus hemolyticus
Purity	> 95% pure
Conjugate	Unconjugated
Size	1 mg
Buffer	Supplied in liquid form dissolved in PBS pH 7.4
Preservative	None
Storage	Aliquot and store at -20 or -80°C. Avoid repeated freeze/thaw cycles.

BACKGROUND

Introduction	Streptolysin is a streptococcal hemolytic exotoxin. Types include Streptolysin O (SLO), which is oxygen-labile, and Streptolysin S (SLS), which is oxygen-stable.
Keywords	S. hemolyticus Streptolysin O; Streptolysin O, Cal.; Streptolysin O, Calibrator; Streptolysin O; Thiol activated cytolysin; Slo; S. hemolyticus Streptolysin O; Streptococcus hemolyticus Streptolysin O; S. hemolyticus; Streptococcus hemolyticus