



3-Methyl-Quinoxaline-2-Carboxylicacid [HRP] (DAGA-044H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	3-Methylquinoxaline-2-carboxylic Acid;3-Methyl-Quinoxaline-2-Carboxylicacid
----------	---