



Clonazepam [KLH] (DAGA-098K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------|---|
| Conjugate | KLH |
| Applications | Immunogen |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|----------|---|
| Keywords | Clonazepam;Rivotril;Anteplepsin;Iktorivil |
|----------|---|