



## L. pneumophila Antigen (DAGA-165)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Purified native Legionella pneumophila protein
<b>Antigen Description</b>	Legionella pneumophila is a thin, aerobic, pleomorphic, flagellated, non-spore forming, Gram-negative bacterium of the genus Legionella. L. pneumophila is the primary human pathogenic bacterium in this group and is the causative agent of Legionnaires' disease, also known as legionellosis.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB standard, antibody ELISA, immunogen, etc.
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	50 mM glycine buffer, pH 9.0
<b>Preservative</b>	None
<b>Storage</b>	≤-70°C

### BACKGROUND

**Keywords** L. pneumophila; Legionella pneumophila; L. pneumophila Oligoribonuclease; L. pneumophila orn; Legionella pneumophila Oligoribonuclease; Legionella pneumophila orn; L. pneumophila (pal);