



Creatinine [BTG] (DAGA-186G)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Creatinine is a breakdown product of creatine phosphate in muscle, and is usually produced at a fairly constant rate by the body.
Conjugate	BTG
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.5 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Creatinine Acid; Creatine phosphate
----------	-------------------------------------