



RRV Envelope Antigen 2 (aa 271-334) [His] (DAGA-262)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Ross River virus belongs to the virus genus Alphavirus, which is part of the family Togaviridae. The alphaviruses are a group of small enveloped single-strand positive-sense RNA viruses. |
| Purity | >95% , based on SDS PAGE |
| Conjugate | His |
| Applications | WB standard, antibody ELISA, immunogen, etc. |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 100 µg, 1 mg |
| Buffer | In PBS with 8M Urea |
| Preservative | None |
| Storage | Stable for 6 months at -20°C from the date of shipment. Nonhazardous. No MSDS required. |

BACKGROUND

| | |
|-----------------|---|
| Keywords | RRV E2;Ross River Virus envelope protein;RRV Envelope Antigen 2;RRV |
|-----------------|---|