



RRV Nonstructural 2 (aa 535-1332) [His, GST] (DAGA-263)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Ross River virus belongs to the virus genus Alphavirus, which is part of the family Togaviridae. The alphaviruses are a group of small enveloped single-strand positive-sense RNA viruses.
Purity	>90% , based on SDS PAGE
Conjugate	His, GST
Applications	WB standard, antibody ELISA, immunogen, etc.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg, 1 mg
Buffer	In PBS with 8M Urea
Preservative	None
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.

BACKGROUND

Keywords RRV nsp2 GST;Ross River Virus non-structural protein 2;RRV nsp2;RRV Nonstructural 2;RRV