



RRV Nonstructural 3 (aa 1333-1875) [His, GST] (DAGA-264)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Ross River virus belongs to the virus genus Alphavirus, which is part of the family Togaviridae. The alphaviruses are a group of small enveloped single-strand positive-sense RNA viruses. |
| Purity | >90% , based on SDS PAGE |
| Conjugate | His, GST |
| Applications | WB standard, antibody ELISA, immunogen, etc. |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 100 µg, 1 mg |
| Buffer | In PBS with 8M Urea |
| Preservative | None |
| Storage | Stable for 6 months at -20°C from the date of shipment. Nonhazardous. No MSDS required. |

BACKGROUND

Keywords RRV nsp3 GST;Ross River Virus non-structural protein 3;RRV nsp3;RRV Nonstructural 3;RRV