



Mayaro virus E1 glycoprotein (aa 1-413) [His] (DAGA-270)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Mayaro Virus has a structure similar to other Alphavirus. It is an enveloped virus and has an icosahedral capsid. Infection with Mayaro virus causes an acute, self-limited dengue-like illness of 3–5 days' duration.
Purity	>95% , based on SDS PAGE
Conjugate	His
Applications	WB standard, antibody ELISA, immunogen, etc.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg, 1 mg
Buffer	In PBS with 8M Urea
Preservative	None
Storage	Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Nonhazardous. No MSDS required.

BACKGROUND

Keywords	MV E1;Mayaro virus E1;Mayaro Virus envelope glycoprotein 1;Mayaro virus E1 glycoprotein;Mayaro virus;MAYV
-----------------	-----------------------------------------------------------------------------------------------------------