



Native Streptococcus Pneumoniae (DAGA-3026)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Suitable for use as a positive control in Lateral Flow assays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Lysate
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.1 ml
Buffer	137 mM Sodium Chloride (NaCl), 2.7 mM Potassium Chloride (KCl), 10 mM Disodium Phosphate (Na ₂ HPO ₄), 1.8 mM Dihydrogen Phosphate (H ₂ PO ₄), pH 7.4 + 0.2.
Preservative	0.09% Sodium Azide
Storage	Store for less than 1 week at 2-8°C. For long term storage, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.