



## Native Scl-70 Antigen (DAGA-3028)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Conjugate</b>	N/A
<b>Applications</b>	Standard solid phase immunoassays. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1000 UNITS
<b>Buffer</b>	0.01 M Tris, 0.53 M Sodium Chloride.
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at $-70^{\circ}\text{C} + 10^{\circ}\text{C}$ . Aliquot to avoid multiple freeze/thaw cycles.

### BACKGROUND

<b>Keywords</b>	Scl-70; DNA topoisomerase I; Topo I; TOP1; Systemic sclerosis
-----------------	---