



Native West Nile Virus antigen (Strain 385-99 (New York)) (DAGA-3052)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	PBS buffer, pH 7.4 with <0.1% CHAPS
Preservative	None
Storage	Store at -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	West nile virus; WNV
-----------------	----------------------