



## Native West Nile Virus antigen (Strain B956 (Uganda)) (DAGA-3053)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Conjugate</b>	N/A
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	PBS buffer, pH 7.4 with <0.1% CHAPS
<b>Preservative</b>	None
<b>Storage</b>	Store at -80°C. Avoid multiple freeze/thaw cycles.

### BACKGROUND

<b>Keywords</b>	West nile virus; WNV
-----------------	----------------------