



Native Respiratory Syncytial Virus (Strain Long) Antigen (DAGA-3443)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	ELISA
Format	Liquid
Concentration	1 mg/mL
Size	1 mg
Buffer	50 mM Glycine, 140 mM NaCl, pH 9.0
Storage	Store at -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	Respiratory syncytial virus; RSV; RSV A
----------	---