



Recombinant phospholipase B-like 2 protein (PLBL-2) [His] (DAGF-124)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant PLBL-2 (UniProt #G3I6T1) with 6xHIS tag was expressed in and purified from CHO cells. Due to the processing mechanism of PLBL-2, the protein is detected at estimated molecular weights of 70 kDa and 44 kDa under reducing conditions.
Purity	The product is determined to be >85% pure based on HPLC analysis.
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg
Buffer	PBS buffer, pH 7.2
Preservative	None
Storage	Store at -80°C. After initial thaw, prepare multiple aliquots and store at -80°C. Multiple freeze-thaw cycles are not recommended.

BACKGROUND

Keywords phospholipase B-like 2 protein; PLBL-2

GENE INFORMATION

UniProt ID UniProt #G3I6T1