



Recombinant HIV type 1 TAT Protein [FITC] (DAGH003)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	E.coli derived tat (HIV type 1) conjugated to FITC
Specificity	This fluoresce conjugated tat binds to murine monoclonal antibodies of defined epitope specificity.
Conjugate	FITC
Applications	Dot Blots; FACS; Immunohistocytology; Direct Immunofluorescence
Stability	At least 6 months at -75°C
Bio-activity	Not determined
Format	Frozen Liquid
Size	50 µg
Preservative	None
Storage	-75 °C

BACKGROUND

Keywords HIV type 1 tat; HIV-1 tat; human immunodeficiency virus-1 tat; human immunodeficiency virus; 30 kDa HIV-1 TAT-interacting protein; HIV-1 TAT-interactive protein 2; Oxidoreductase HTATIP2